


CAT. NO.....	296-3345	296-4545
PORE SIZE (um).....	0.45	0.45
UPPER CAP. (mL).....	150	150
BOTTLE NECK SIZE.....	33mm	45mm
MEMBRANE COLOR.....	WHITE	WHITE
MEMBRANE DIA (mm)...	50	50
	INCHES	MM
A (UPPER I.D.).....	3.15	80.0
B HEIGHT.....	3.6	91
C DIA (BASE).....	1.4	36
D DIM.....	2.5	64
NOM. WT. (gms).....	38	40

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	BLUE PS
ADAPTOR	WHITE LOPE
GASKET	FOAMED LOPE
MEMBRANE	PES
VENT PLUG	CA

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Pancromo Creek Drive
Rochester, NY 14625-2385, USA

NAME	MC ²	PVC	PSS 00-0296-XX45	REV.	1205
DR'N	SH	12/22/05			
DESCRIPTION:					

BOTTLE TOP FILTERS

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.