


CAT. NO.....	295-3345	295-4545
PORE SIZE (um).....	0.45	0.45
UPPER CAP. (mL).....	500	500
BOTTLE NECK SIZE.....	33mm	45mm
MEMBRANE COLOR.....	WHITE	WHITE
MEMBRANE DIA (mm)...	75	75
	INCHES	MM
A (UPPER I.D.).....	3.8	96
B HEIGHT.....	5.0	127
C DIA (BASE).....	1.4	36
D DIM.....	2.5	64
NOM. WT. (gms).....	70	72

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	BLUE PS
ADAPTOR	WHITE LOPE
GASKET	FOAMED LOPE
MEMBRANE	PES
VENT PLUG	CA

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Pancrome Creek Drive
Rochester, NY 14625-2385, USA

MFG MC² PVC
DR'N SH 12/22/05
DESCRIPTION:

PSS 00-0295-XX45

REV 1205

BOTTLE TOP FILTERS

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.