


CAT. NO.	290-3320	290-3345	290-4520	290-4545
PORE SIZE (µm)	0.2	0.45	0.2	0.45
UPPER CAP. (mL)	150	150	150	150
BOTTLE NECK SIZE	33mm	33mm	45mm	45mm
MEMBRANE COLOR	WHITE	WHITE	WHITE	WHITE
MEMBRANE DIA (mm)	50	50	50	50
	INCHES	MM	INCHES	MM
A (UPPER I.D.)	3.15	80.0	3.15	80.0
B HEIGHT	3.6	91	3.6	91
C DIA (BASE)	1.4	36	1.9	48
D DIM.	2.5	64	3.0	76
NOM. WT. (gms)	38	38	40	40

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	YELLOW PS
ADAPTOR	WHITE LOPE
GASKET	FOAMED LOPE
MEMBRANE	SFCA
VENT PLUG	CA

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Pancromo Creek Drive  
Rochester, NY 14625-2385, USA

NAME	MC <sup>2</sup> PVC	PSS 00-0290-XXXX	REV	0602
DATE	8/27/02			
DESCRIPTION:	BOTTLE TOP FILTERS			

CONFIDENTIAL  
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.