


CAT. NO.	292-3320	292-4520	297-3320	297-4520
PORE SIZE (µm)	0.2	0.2	0.2	0.2
UPPER CAP. (mL)	1000	1000	1000	1000
BOTTLE NECK SIZE	33mm	45mm	33mm	45mm
MEMBRANE COLOR	WHITE	WHITE	WHITE	WHITE
MEMBRANE DIA (mm)	90	90	90	90
	INCHES	MM	INCHES	MM
A (UPPER I.D.)	4.49	114.0	4.49	114.0
B HEIGHT	6.43	163.3	6.43	163.3
C DIA (BASE)	1.39	35.3	1.39	35.3
D DIM.	2.48	63.0	2.48	63.0
NDM. WT. (gms)	117	118	117	118

MATERIALS	
CATALOG NUMBER	292-3320/-4520 297-3320/-4520
UPPER CONTAINER	PS
COVER	PS
COLLAR	YELLOW PS
ADAPTOR	WHITE LOPE
GASKET	FOAMED LOPE
MEMBRANE	SFCA
VENT PLUG	CA

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Pancromo Creek Drive  
Rochester, NY 14625-2385, USA

NAME: MC<sup>2</sup>  
CR: B  
DESCRIPTION:

APP: 0

PSS 00-029X-XX20

REV: 0902

BOTTLE TOP FILTERS

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.