

TECHNICAL DATA SHEET

Part Number	028649
Description	Unwired rack for tubes



Benefits :

- ✓ Securely hold a full load of test tubes or centrifuge tubes for 30 mm Test Tube Racks.
- ✓ Molded in a single piece for exceptional strength and durability
- ✓ These racks have no coating to chip, peel or rust

Certificate :

- ✓ Autoclavable

Technical specifications :

- ✓ For use with 30 mm tubes
- ✓ Color : white
- ✓ Shape : rectangle
- ✓ Shelves : 3
- ✓ Diameter inner : 30 mm
- ✓ Dimensions : 28,3 x 10,8 x 9,2 cm

Sales unit :

1 Unwired rack for tubes