


CATALOG NO.	NOMINAL TUBE SIZE, (ml)	MATERIAL
005973	0.5	RESMER

TUBE SIZE.....	0.5	1.5
ARRAY.....	4 x 6	4 x 6
	INCHES	INCHES
A (LENGTH).....	5.73±.16	145.5±4.0
B (WIDTH).....	3.62±.14	91.9±3.6
C (HEIGHT).....	2.15±.10	54.6±2.5
D (DIA.).....	0.33±.01	8.4±0.3
NDM. WT. (gn).....	85	11.2±0.3
COLOR	SIZE CODE	
WHITE	-0005	-0015

TUBE SIZE.....	0.5	1.5
ARRAY.....	8 x 12	8 x 12
	INCHES	INCHES
A (LENGTH).....	6.75±.16	171.5±4.0
B (WIDTH).....	4.42±.14	112.3±3.6
C (HEIGHT).....	2.15±.10	54.6±2.5
D (DIA.).....	0.33±.01	8.4±0.3
NDM. WT. (gn).....	105	140
COLOR	SIZE CODE	
WHITE	-9005	-9015
BLUE	-9305	-9315
RED	—	-9515

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Poncaque Creek Drive
Rochester, NY 14625-2385, USA

MFG	MC	APP	REV
OK'U	SHARAD	1/21/03	0103

PSS 00-5973-XXXX-1

MICROCENTRIFUGE TUBE RACKS

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.