


CATALOG NO.	NOMINAL CAPACITY, (mL)		MATERIALS			
	TUBE	GASKET	CLOSURE	TUBE	GASKET	CLOSURE
003137	—	30	50	PSF	RED SILICONE	WHITE PP
003138	10	16	30	PC	RED SILICONE	WHITE PP

SIZE CODE.....	-0010	-0016		
NOM. CAP. (mL).....	10	16		
CLOSURE SIZE (mm).....	13	15		
NOM. WEIGHT W/ GASKET & CLOSURE.....	6	10		
	INCH	MM	INCH	MM
A (NECK I.O.).....	0.32	8.1	0.38	9.6
B (TUBE O.D.).....	.635 ⁺⁰ _{-.02}	16.13 ⁺⁰ _{-.05}	.715 ⁺⁰ _{-.03}	18.16 ⁺⁰ _{-.076}
C (HT. W/O CLOSURE).....	3.12	79.2	3.86	98.0
D (HT. W/ CLOSURE).....	3.30	83.8	4.22	107.2

SIZE CODE.....	-0030	-0050		
NOM. CAP. (mL).....	30	50 *		
CLOSURE SIZE (mm).....	20	24		
NOM. WEIGHT W/ GASKET & CLOSURE.....	14	24		
	INCH	MM	INCH	MM
A (NECK I.O.).....	0.54	13.7	0.64	16.3
B (TUBE O.D.).....	1.016 ⁺⁰ _{-.04}	25.81 ⁺⁰ _{-.102}	1.135 ⁺⁰ _{-.028}	28.83 ⁺⁰ _{-.064}
C (HT. W/O CLOSURE).....	3.63	92.2	4.09	103.9
D (HT. W/ CLOSURE).....	4.04	102.6	4.52	114.8

* ACTUAL BRIM CAPACITY OF -0050 IS 43mL

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CAT. NO. 3138-0010

CAT. NO. 3138-0016
3137-0030
3138-0030
3137-0050
3138-0050
(TYP)

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MAKE	TL	APP'D	REV.
CF-0	OV	GRP	1009

DESCRIPTION:
PSS 00-313X-00XX-1

DAK RIDGE CENTRIFUGE TUBES
WITH SEALING CAP