


CATALOG NUMBER.....	6379-0010	
NOM. CONTAINER CAPACITY (L).....	10	
CLOSURE SIZE (mm).....	83B	
NOM. WT. (gms).....	1920	
	INCH	MM
FUNNEL TOP I.D.....	10	254
NOM. OVERALL HT. W/ LID CLOSED.....	20.5	521
NOM. OVERALL WIDTH W/ LID CLOSED.....	11.4	290

COMPONENT	MATERIAL
COVER	RED HDPE
FUNNEL	HDPE
NECK SCREEN INSERT	PP
83B CLOSURE	WHITE PP
FUNNEL ADAPTER	PP
FILTER HOUSING/MEMBRANE	PP/PTEE
10L CARBOY W/ CLOSURE	FLPE/WHITE PP

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
S C I E N T I F I C

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

NAME	MCC	PVC	PVC
DATE	PVC	11/30/01	
DESCRIPTION:	PSS 00-6379-0010		
REV.	1101		

10L SAFETY WASTE FUNNEL SYSTEM