


CATALOG NO.	4301-2000	
SIZE (NOM VOL. mL)	2000	
NOM. BRIM CAP. (mL)	2400	
NOM. WEIGHT (gms)	280	
	INCH	MM
A (NECK I.D.)	1.53 ± 0.02	38.9 ± 0.5
B (HEIGHT W/ CLOSURE)	15.8 ± 0.2	403.9 ± 5.1
C (HEIGHT W/O CLOSURE)	15.6 ± 0.2	396.2 ± 5.1
D (O.D. AT BULB AREA)	6.21 ± 0.09	157.7 ± 2.0
E (TUBING O.D.)	0.34 ± 0.01	8.7 ± 0.3
W.T. (NOM. WALL THK.)	0.04	1.0

MATERIALS	
CLOSURE	TEFZEL ETFE
FUNNEL BODY	TEFLON FEP
STOPCOCK ASSEMBLY	TEFLON TFE
TIP	TEFLON FEP

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Pancrome Creek Drive  
Rochester, NY 14625-2385, USA

NAME	MFC	APP'D	PVC	PSS 00-4301-2000-S	REV	0592
DATE	PVC	5-4-92				

DESCRIPTION:

SEPARATORY FUNNEL

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.