


CATALOG NO.	5245-0040	5245-0050
SIZE CODE	D	E
MAX. PIPET LENGTH (MM/IN.)	406/16	610/24
BASE	BLUE HOPE	
BODY	HOPE	
TUBING *	PVC	
	INCHES	MM
A (BASE DIA.)	9.0	229
B (OVERALL WIDTH)	10.7	272
D (NOM. OVERALL HEIGHT)	30	762
H (NOM. HEIGHT)	27	686
NOM. WT. (gms)	1800	1700

* INCLUDES (2) PIECES OF TUBING AS FOLLOWS:
 (1) 1" I.D. x 1 1/4" O.D. x 18" LONG
 (1) 1/2" I.D. x 3/4" O.D. x 48" LONG

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
 75 Pancromo Creek Drive
 Rochester, NY 14625-2385, USA

MFG	MCC	APP'D	PVC
OK'D	PVC	9/14/95	
DESCRIPTION:			

PSS 00-5245-00X0

REV 0995

PIPET WASHER-RINSER

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.