


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS			
		BOTTLE	CLOSURE	PLUG	O-RING
3140	1000	PC	PPE+HIPS	PP-GLASS FILLED	SILICONE

SIZE CODE.....	1002/1006/1024	
NOM. CAP.....	1000	
CLOSURE SIZE (mm) .....	84.5	
	INCH	MM
A (NECK I.D.).....	2.82±.02	71.6±0.5
B (HT. W/ CLOSURE).....	7.69±.06	195.3±1.5
C (DIA.).....	3.84±.05	97.5±1.3
NOM WALL THICKNESS:..	.08	2.0

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

**ThermoFisher**  
S C I E N T I F I C

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

NAME: PB GRP  
CUST: KOW 8/24/07  
DESCRIPTION:

PSS 003140-10XX

REV: 0807

1L HIGH SPEED CENTRIFUGE BOTTLE

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET  
CONTAINS CONFIDENTIAL AND/OR  
PROPRIETARY INFORMATION AND  
IS TO BE USED ONLY FOR THE  
EXPRESSED PURPOSES AGREED UPON  
BY AN AUTHORIZED OFFICER OF  
THERMO FISHER SCIENTIFIC.