

ratiolab® Waste Disposal Bag

ratiolab® Waste Disposal Bag, Standard

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 134 °C, for autoclaving, the bag must not be completely closed
- only for bag size 200 x 300 mm: practical bench-top stand, plastic-coated, with 3 rubber feet



W x H mm	Volume approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
200 x 300	2.5	50 µm	dispenser 1 x 100	100	70 01 000
300 x 500	8	50 µm	cardboard box 1 x 500	500	70 01 100
400 x 780	38	50 µm	cardboard box 1 x 500	500	70 01 150
600 x 800	60	50 µm	cardboard box 1 x 500	500	70 01 200
600 x 800	60	50 µm	cardboard box 1 x 100	100	70 01 210
600 x 800	60	100 µm	cardboard box 1 x 125	125	70 01 212
700 x 1100	110	50 µm	cardboard box 1 x 350	350	70 01 300
700 x 1100	110	50 µm	cardboard box 1 x 75	75	70 01 310
700 x 1100	110	100 µm	cardboard box 1 x 75	75	70 01 311
bench-top stand for waste disposal bags 200 x 300 mm				1	70 01 400
Set consisting of bench-top stand and 100 waste disposal bags (standard) 200 x 300 mm				1	70 01 401

ratiolab® Waste Disposal Bag, High Transparency

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 121 °C, for autoclaving, the bag must not be completely closed
- only for bag size 200 x 300 mm: practical bench-top stand, plastic-coated, with 3 rubber feet



W x H mm	Volumen approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
200 x 300	2.5	50 µm	dispenser 1 x 100	100	70 01 005
300 x 500	8	50 µm	cardboard box 1 x 500	500	70 01 105
400 x 780	38	50 µm	cardboard box 1 x 500	500	70 01 155
600 x 800	60	50 µm	cardboard box 1 x 500	500	70 01 205
700 x 1100	110	50 µm	cardboard box 1 x 350	350	70 01 305
bench-top stand for waste disposal bags 200 x 300 mm				1	70 01 400
Set consisting of bench-top stand and 100 waste disposal bags (high transparency) 200 x 300 mm				1	70 01 405