


## 20 $\mu$ l (1-20 $\mu$ l) Barrier - Catalog #30370T



### Catalog #

30370T

### Description

Natural, Polypropylene, Low Binding, Sterile

### Packaging Type

Racked Tips

### Packaging Details

Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, Hinged Rack

### Attributes



Certified RNase/DNase Free  
Certified Human DNA Free  
Certified PCR Inhibitor Free  
Certified Pyrogen Free

- **Material:** Polypropylene
- **Volume Range:** 1-20  $\mu$ l
- **Grad Lines:** No Graduations