


200 µl (1-200 µl) Standardization - Catalog #30430



Catalog #	Description
30430	Natural, Polypropylene
Packaging Type	Packaging Details
QuickRack	Non-Sterile, Racked, 960 Tips/Stack, 5 Stacks/Case, QuickRack
Attributes	



Certified Human DNA Free  
 Certified PCR Inhibitor Free  
 Certified Pyrogen Free

**Volume Range:** 1-200 µl  
**Grad Lines:** 10, 50, 100 µl

DOMINIQUE DUTSCHER SAS