

TECHNICAL DATA SHEET

Our flexible polystyrene boats feature smooth surfaces and rounded corners to provide accurate pour-outs with minimal sample loss. Boats have a diamond shape that easily flexes for accurate pouring. Feature a flat bottom to resist tipping. Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases. May be used with microwave moisture analyzers. Pack of 500.

Product name:	Diamond-Shaped Weighing Boats
Part number:	HS1426B
Type:	Black
Material of construction:	Polystyrene
Dimensions:	78 x 56 x 14 mm
Volume:	30 mL
Resist temperatures up to:	158°F (70°C)
Antistatic	Yes