



TECHNICAL DATA SHEET

Our flexible polystyrene boats feature smooth surfaces and rounded corners to provide accurate pour-outs with minimal sample loss. Boats are shallow with a wide, flat bottom to resist tipping. Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases. May be used with microwave moisture analyzers. Pack of 500.

Product name:	Standard Weighing Boats
Part number:	HS1423A
Type:	Black
Antistatic	Yes
Material of construction:	Polystyrene
Material additive:	
Dimensions:	46 x 46 x 8 mm
Size:	Small
Resist temperatures up to:	158°F (70°C)

Heathrow Scientific®, LLC

620 Lakeview Parkway • Vernon Hills • Illinois • 60061 • Phone: (847) 816-5070 • Fax: (847) 816-5072
E-mail: info@heathrowscientific.com • Website: www.heathrowscientific.com