


ATC

Separate temperature sensor

°C ISM



Sensor type temperature probe

Order number 51300164

Specifications

pH range

zero point (pH) -

temperature range -30...100 °C

temperature probe NTC 30 KΩ

shaft material stainless steel

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.0 m

connector Cinch

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

Dry