


InLab[®] 731-ISM

Cond °C ISM



For general applications like standard and high conductivity samples.



Sensor type conductivity cell

Order number 30014092

Specifications

measuring range 0.01...1000 mS/cm

zero point (pH) -

temperature range 0...100 °C

temperature probe NTC 30 kΩ

shaft material epoxy

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.2 m

connector Mini-DIN

cell constant (cm⁻¹) 0.57

cell type 4 graphite poles

indicating metal -

Storage

dry