


InLab® Solids Pro-ISM

For solids and semisolids samples like meat, cheese, fruits, vegetables

pH Ref °C ISM



Scale 1:1



Sensor type combined pH electrode

Order number 51344155

Specifications

pH range 1...11
zero point (pH) 7.00 ± 0.25

temperature range 0...80 °C
temperature probe -

shaft material glass
membrane glass Lot
membrane resistance < 250 MΩ

reference system ARGENTHAL™
reference diaphragm open junction
reference electrolyte XEROLYT® EXTRA polymer
bridge electrolyte -

electrode head MultiPin™
cable length -
connector -

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage

InLab storage solution (30111142)
or FRISCOLYT-B® (51340053)