


# InLab<sup>®</sup> Viscous

For viscous samples like cosmetics, paints, glue

pH Ref °C ISM



Scale 1:1



**Sensor type** combined pH electrode

**Order number** 51343150

**Specifications**

pH range 0...14

zero point (pH) 7.00 ± 0.25

temperature range 0...130 °C

temperature probe -

shaft material glass

membrane glass HA

membrane resistance < 600 MΩ

reference system SteadyForce™

reference diaphragm ceramic

reference electrolyte FRISCOLYT-C®

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm<sup>-1</sup>) -

cell type -

indicating metal -

**Storage**

InLab storage solution (30111142)

or 3 mol/L KCl (51350072)