


LE427 IP67

For semi-solid and solid samples as
meat, fruits, cheese and soil

pH Ref °C ISM



Sensor type combined pH electrode

Order number 30259840

Specifications

pH range 1...11
zero point (pH) 7.00 ± 0.25

temperature range 0...50 °C
temperature probe -

shaft material glass
membrane glass LoT
membrane resistance < 250MΩ

reference system Ag / AgCl
reference diaphragm ceramic
reference electrolyte gel
bridge electrolyte -

electrode head -
cable length 2.0 m
connector BNC

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage

3M KCl (51350072)