


InLab[®] 741

Cond °C ISM



For low conductivity range samples
like pure water and ultra pure water.



Sensor type conductivity cell

Order number 51344024

Specifications

measuring range 0.001...500 $\mu\text{S}/\text{cm}$
zero point (pH) -

temperature range 0...100 °C
temperature probe NTC 30 k Ω

shaft material stainless steel
membrane glass -
membrane resistance -

reference system -
reference diaphragm -
reference electrolyte -
bridge electrolyte -

electrode head -
cable length 1.2 m
connector Mini-DIN

cell constant (cm^{-1}) 0.105 (certified)

cell type 2 steel poles
indicating metal -

Storage

dry