

## TECHNICAL DATA SHEET

---

<b>Part Number</b>	033283
<b>Description</b>	Microtube



### Benefits :

- ✓ Made of virgin polypropylene
- ✓ Manufactured without mold release agents, plasticizers or biocides: avoids extractables and guarantees the integrity of the samples
- ✓ Safe-Lock: locking system
- ✓ Centrifugeable up to 30,000 g in angle rotor
- ✓ Operating temperatures: from -80°C to 100°C
- ✓ Autoclavable (121°C, 20 min)
- ✓ Natural color
- ✓ Volume 0.5ml

### Sales unit :

500 Microtubes

---