

TECHNICAL DATA SHEET

Part Number	033520
Description	Microtube PCR 0,5 ml



Benefits :

- ✓ Thin-walled Eppendorf microtube specially designed for PCR applications
- ✓ Colours: natural, high transparency
- ✓ With marking surface
- ✓ Without graduations
- ✓ With hinged cap, easy opening Excellent sealing
- ✓ Non-sterile, free of DNases, RNases
- ✓ PCR-clean

Technical specifications :

- ✓ For use with 0.5 ml thermoblocks
- ✓ Centrifuge up to 10,000 g
- ✓ Volume: 0.5ml

Sales unit :

500 Microtubes PCR
