

TECHNICAL DATASHEET

PART NUMBER

035100ACL



The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.

**RNase
DNase**

**Endo-
toxins**


**Human
DNA**

Pyrogen

ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	035100ACL
Packaging	ClearLine rack with slip lid
Volume	200 µl
Graduation	10, 50, 100
Length	51,44 mm
Color	Natural
Filter	No
Low Retention	Yes
Sterility	Yes
Autoclavable	Yes



Sales units :

960 tips