


20 μ l (1-20 μ l) Barrier - Catalog #35220



Catalog

35220

Description

Natural, Polypropylene, Low Binding, Sterile

Packaging Type

Racked Tips

Packaging Details

Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, Cell Rack

Attributes



Certified RNase/DNase Free
Certified Human DNA Free
Certified PCR Inhibitor Free
Certified Pyrogen Free

- **Material:** Polypropylene
- **Volume Range:** 1-20 μ l
- **Grad Lines:** No Graduations