

## EP-2160 | RUNONE™ SYSTEM, TIMER AND MULTICASTER, ORANGE, 200-240V

includes: [EP-2002](#), [EP-1139](#), [EP-1010](#), [EP-1001](#), [EP-1017](#)



### **RunOne Tank, Qty:1**

Catalog Number: EP-2002

Material: White Polycarbonate

### **RunOne Power Supply, Qty:1**

Catalog Number: EP-1139, RunOne Power Supply with Timer

Material: White Polycarbonate

Dimensions: LxWxH: 19 cm x 12.5 cm x 5.5 cm (7.5" x 5" x 2.25")

Maximum Voltage Input: 120VAC, 60Hz

Over Current Limit: 300mA

Power Rating: 70W

Electrical Rating: 120VAC, 60Hz, 70W, 600mA

### **Transformer, Qty:1**

Catalog Number: EP-1010

Voltage Input: 230V~ 50Hz 60VA

Voltage Output: 120V~ 50Hz 50VA

### **RunOne Casting System, Qty:1**

Catalog Number: EP-1001

Casting Stand Material: White Polycarbonate

Tray Material: Clear Polycarbonate w/ Black Print

Mini Trays: 5.4(w) x 5.9(l)cm, Qty:4

Landscape Trays: 10.9(w) x 5.9(l)cm, Qty:2

Combs - 1mm thick: 6 and 8 wells for Mini; 12 and 17 wells for Landscape; Qty:2

Load Volume: 28µl for 6 and 12 well, 17µl for 8 and 17 wells

### **MultiCaster System, Qty:1**

Catalog Number: EP-1017, Orange MultiCaster System

Material: Orange Polycarbonate

Casting Stand: Orange Casting Stand, 2 Slots, Qty:1

Gel Tray: Med UV Transparent 11.7(w) x 8.9(l)cm (EP-1024), Qty:1

Combs: Orange 2+1 (250µl)/5+2 (3x 75µl, 2x 112µl), Qty:2

MultiChannel Compatible: No

Visualization Plate: Red Med Plate 11.5(w) x 8.9(l)cm (EP-1048), Qty:1

Gel Size: 11.7(w) x 8.9(l) cm