

# Cryo-Gloves®

Cryogenic Protection in Ultra-Cold Environments



10% of PINK sales go to Cancer Research



GLOVE LENGTH	WRIST	MID-ARM	ELBOW	SHOULDER
	11"-13" (280-330mm)	13¼"-15½" (335-395mm)	17¼"-19¾" (440-500mm)	24½"-27¼" (620-695mm)
GLOVE SIZE				
<b>Small/8</b>	Blue ● WRS Pink ● P-WRS	Blue ● MAS Pink ● P-MAS	Blue ● EBS Pink ● P-EBS	Blue ● SHS Pink ● P-SHS
<b>Medium/9</b>	Blue ● WRM Pink ● P-WRM	Blue ● MAM Pink ● P-MAM	Blue ● EBM Pink ● P-EBM	Blue ● SHM Pink ● P-SHM
<b>Large/10</b>	Blue ● WRL Pink ● P-WRL	Blue ● MAL Pink ● P-MAL	Blue ● EBL Pink ● P-EBL	Blue ● SHL Pink ● P-SHL
<b>X-Large/11</b>	Blue ● WRXL Pink ● P-WRXL	Blue ● MAXL Pink ● P-MAXL	Blue ● EBXL Pink ● P-EBXL	Blue ● SHXL Pink ● P-SHXL
<b>XX-Large/12</b>	Blue ● WRXXL Pink ● P-WRXXL	Blue ● MAXXL Pink ● P-MAXXL	Blue ● EBXXL Pink ● P-EBXXL	Blue ● SHXXL Pink ● P-SHXXL
<b>3X-Large/13</b>	Blue ● WR3XL Pink ● P-WR3XL	Blue ● MA3XL Pink ● P-MA3XL	Blue ● EB3XL Pink ● P-EB3XL	Blue ● SH3XL Pink ● P-SH3XL

## Features

- Cryogenic protection for ultra-cold applications down to -196°C (-320°F)
- State-of-the art materials and a multi-layered construction allow for a maximum level of thermal protection, flexibility, and dexterity which are essential features when function is important and safety is critical
- High performance, thermal inner lining wicks moisture away from hands, maintaining comfort over extended periods
- Safety certified to meet EU standards: EN 511, EN 388, EN 420
- **Applications:** handling dry ice, low and ultra-low temperature freezers, closed cryogenic systems



CE Cat. III  
Pink Glove C E Pending

**CAUTION** Not for immersion in liquid cryogenics