

# Cryo-Apron®

Protection from Splashes and Contact Cold

## Features

- Cryogenic protection for ultra-cold applications down to -196°C (-320°F)
- Made from state-of-the-art materials that are lightweight for comfort and ease of movement
- Multiple lengths available for maximum safety
- Adjustable straps at the neck and waist to ensure a proper fit
- Quick-release buckles at the neck and waist for swift on/off function
- Safety certified to meet EU standards: EN 511, EN 388, EN ISO 13688:2013
- **Applications:** dispensing/transferring cryogenic liquids, handling of cryogenic samples into or removing from storage vessels



| APRON LENGTH | 36"<br>(91.4cm)              | 42"<br>(106.7cm)             | 48"<br>(121.9cm)             | 54"<br>(137.2cm)             |
|--------------|------------------------------|------------------------------|------------------------------|------------------------------|
| Apron        | Blue ● CA36<br>Pink ● P-CA36 | Blue ● CA42<br>Pink ● P-CA42 | Blue ● CA48<br>Pink ● P-CA48 | Blue ● CA54<br>Pink ● P-CA54 |



CE Cat. III  
Pink Apron C E Pending