

TECHNICAL DATA SHEET

Part Number	035970
Description	Eppendorf SmartBlock 5 mL



Benefits :

- ✓ Thermoblock for 8 Eppendorf Tubes® 5.0 mL
- ✓ Optimized block design supports optimal temperature transfer within your sample to receive excellent results
- ✓ Eppendorf SmartBlock-specific calibration by integrated sensor technology maximizes temperature accuracy and homogeneity
- ✓ Every SmartBlock is individually inspected in production for temperature accuracy and is delivered incl. s serial-number specific certificate for documentation purpose
- ✓ Insulated thermoblocks for ergonomic operation - you won't burn your fingers (PhysioCare Concept)
- ✓ Automatic block recognition for safe handling
- ✓ Eppendorf QuickRelease technology enables fast and simple block exchange - no need for tools

Sales unit :

1 Thermoblock