


100 µl (1-100 µl) Barrier - Catalog #36060



Catalog #	Description
36060	Natural, Polypropylene, Sterile

Packaging Type	Packaging Details
Racked Tips	Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, Cell Rack

Attributes



Certified RNase/DNase Free
 Certified Human DNA Free
 Certified PCR Inhibitor Free
 Certified Pyrogen Free

Material: Polypropylene
Volume Range: 1-100 µl
Grad Lines: No Graduations

DOWNLOAD TECHNIQUE DUTSCHER SAS