




100 µl (1-100 µl) Barrier - Catalog #36060T



Catalog #	Description
36060T	Natural, Polypropylene, Low Binding, Sterile
Packaging Type	Packaging Details
Bulk Tips	Sterile, 96 Tips/Rack, 10 Racks/Pack, 5 Packs/Case, Cell Rack
Attributes	



Certified RNase/DNase Free
 Certified Human DNA Free
 Certified PCR Inhibitor Free
 Certified Pyrogen Free

Material: Polypropylene

Volume Range: 1-100 µl

Grad Lines: No Graduations