



TECHNICAL DATA SHEET

Constructed of durable polypropylene, slide boxes hold either 25 or 100 standard slides (75 x 25 mm and 3 x 1"). Hinged lids are held securely closed by latch lock. Stacking feet maximize stability No liner. The 100-place box measures 210 x 169 x 37 mm

Product name:	Slide box, 100 Place, Blue
Part number:	HS15991A
Color:	Blue
Material of construction:	Polypropylene (PP)
Capacity:	100 Place
Dimensions:	210 x 169 x 37 mm

Heathrow Scientific® , LLC

620 Lakeview Parkway • Vernon Hills • Illinois • 60061 • Phone: (847) 816-5070 • Fax: (847) 816-5072
E-mail: info@heathrowscientific.com • Website: www.heathrowscientific.com