

TECHNICAL DATA SHEET

Sturdy polypropylene rack accommodates individual 0.2 mL tubes or strips of 8 or 12 tubes. Tube wells are easily identified with embossed numbers and letters. Rack measures 130 x 98 x 33 mm. Autoclavable. Pack of 5.

Product name:	96-Well PCR [®] Rack
Part number:	HS23461B
Color:	Natural
Wells:	96 wells (5.2 mm wells dia.)
Tubes:	0.2ml
Material of construction:	Polypropylene
Autoclavable:	Yes
Dimensions, module:	130 x 98 x 33 mm
Packaging:	Pack of 5
Packaging Weight:	0.39 kg