



TECHNICAL DATA SHEET

Create as large a rack as you need by linking together multiple racks. Heavy-duty racks feature a unique system of tabs and slots, that facilitate easy connection and sturdy fit. Autoclavable.

Product name:	4-Way Microtube Rack
Part number:	HS29022A
Color:	Natural
Wells:	32 for 0.5 mL tubes, 8.2 mm 32 for 1.5/2.0 tubes, 11.6 mm 12 for 15 mL tubes, 17.5 mm 4 for 50 mL tubes, 30 mm
Tubes:	1.5/2.0 mL, 05 mL, 15 mL, and 50 mL
Material of construction:	Polypropylene
Autoclavable:	Yes
Dimensions, module:	174 x 95 x 52 mm
Packaging:	Pack of 5
Packaging Weight:	0.77 kg





phone 1.847.816.5070
fax 1.847.816.5072
info@heathrowscientific.com

620 Lakeview Parkway
Vernon Hills, IL 60061
www.heathrowscientific.com

DOMINIQUE DUTSCHER S.A.S.

