

TECHNICAL DATASHEET

PART NUMBER

037640TCL


The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.

**RNase
DNase**

**Endo-
toxins**


**Human
DNA**

Pyrogen

ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

Part number	037640TCL
Packaging	ClearLine rack with slip lid
Volume	10 µl
Graduation	2, 5, 10
Length	45,34 mm
Color	Natural
Filter	No
Low Retention	Yes
Sterility	No
Autoclavable	Yes



Sales units :

960 tips