

Certificate of Quality Assurance

We at Heathrow Scientific are committed to developing and manufacturing high-quality, industry-leading products for scientific research. We certify that the following product has been produced and inspected according to our standard quality assurance processes:

**Parafilm[®] Dispenser, AB, Blue
Part Number HS234525B**

Description of Quality Standard: Injection-molded box, lid, and holder of acrylonitrile-butadiene-styrene with non-toxic blue pigment, which are inert or minimally reactive with most common laboratory reagents. Injection-molded cutter of acrylonitrile-butadiene-styrene with non-toxic blue pigment and steel safety razor, which is inert or minimally reactive with most common laboratory reagents. Inspected in the molding facility for cracks, incomplete or distorted shape, or other damage or defects, which would compromise proper function. Individually packaged.

Our mission is to provide conscientiously designed and thoroughly tested products for laboratory professionals in the pure and applied sciences. We strive to maintain the highest quality standards and welcome all inquiries in this regard.

Nicole Smokowicz
Quality Assurance Officer