



Weighted Stand
for bench top use

5
Year
Warranty

Lock Counter Off
to allow use of the marker
without registering more counts

Hand Tally Adapter
allows you to use the eCount
as a digital tally counter

eCount Colony Counter

The eCount is a simple to use colony counter for counting cells in microbiology applications or tissue culture. This unit combines a multi-function electronic counter with a Sharpie® pen for counting and marking to prevent missing or double-counting colonies in a Petri dish.

- US Patent No. D564385 and EP Registered Design No. 705751-0001 featuring a light, ergonomic body that is balanced for easy manipulation
- The unit has convenient easy to use touch controls for multiple counts, sums, counting up or down and turning the unit on/off
- Count up/down verification, enables you to confirm you counts
- Easily verify each count with either an audible beep or LED flash
- Displays totals for up to 32 individual counts or zones
- Quick and easy to change between different colored markers to check different colony counts in the same view
- Multi-functional counter can be turned off allowing the unit be used for writing
- Up to 999 colonies per zone helps simplify the counting process
- Long Life – CR1225 3.0 V Lithium battery



Item No.	L x W x H in	cm	UOM
HS120000	7.4 x 1.2**	17 x 3**	1 ea, 6 ea/cs
HS120228	CR1225 Battery*		1 ea

Package Includes	
eCount Colony Counter	Hand Tally Adapter
Black Sharpie® Pen	Weighted Stand
Long-life Button Battery	

*Battery not included in 5 Year Warranty **Measurements shown are with Sharpie® Pen inserted

