



# TECHNICAL DATA SHEET

Description of Quality Standard: Fabricated of acrylonitrile-butadiene-styrene with no pigment, which is inert or minimally reactive with most common laboratory reagents. Inspected at the fabrication facility for cracks, incomplete or distorted shape, or other damage or defects, which would compromise proper function. Individually packaged.

<b>Product name:</b>	<b>Pipettor Stand, Three-Place</b>
<b>Part number:</b>	<b>HS206203</b>
<b>Color:</b>	<b>Clear</b>
<b>Material of construction</b>	<b>Acrylic</b>
<b>Dimensions:</b>	<b>Slot Upper 25 mm</b> <b>Slot Lower 19 mm</b>
<b>Packaging:</b>	<b>1 ea.</b>

