

Cardboard Slide Folders

- Easily organize slides with the index on the front cover that corresponds with the number grid inside
- Recessed slide compartments and inside covers protect specimens. Slides cannot be touched when the cover is closed, reducing contamination
- Front protection flaps are attached with durable tape for repeated use. Fold flaps flat under the slide holder to save space when in use
- Thumb cut in slide compartments enables easy removal of slides with slight pressure at one end
- Compartment holder is made from a double layer of densely pressed cardboard. Cover is a single layer cardboard
- Available in two sizes that hold 20 or 30 standard (26 x 76mm and 1 x 3in) slides
- Disposable

Item No.	Color	Capacity	L x W x H <i>in</i>	L x W x H <i>cm</i>	UOM
HS9902	Black				
HS9913	Blue				
HS9914	Green	20 slides	8.2 x 13.4 x 0.4	20.8 x 34.1 x 0.9	72/cs
HS120196	Red				
HS120197	Yellow				
HS9930	Black	30 slides	8.2 x 12.9 x 0.4	28.1 x 32.8 x 1	40/cs



DOMESTIC