

TECHNICAL DATA SHEET

Part Number	038338
Description	PCR microtubes 0.2ml



Benefits :

- ✓ Polypropylene
- ✓ Can be used in real-time PCR
- ✓ Colors : natural, "ultra clear" transparency
- ✓ Very thin walls
- ✓ With flat snap cap
- ✓ Autoclavable
- ✓ Certificates: Sterile, free of DNases, RNases, pyrogens and metals E.U. product: extremely uniform
- ✓ "Extra low adhesion" microtube

Technical specifications :

- ✓ Wall thickness: 0.25 mm
- ✓ Evaporation rate: between 2 and 3% (reduced thanks to the homogeneity of the walls)
- ✓ Height: 20.9 – 21.8 mm

Sales unit :

1000 PCR Microtubes
