

TECHNICAL DATA SHEET

Part Number	038366
Description	Thermo Scientific™ ART™ Barrier Extended Length Pipette Tips



Benefits :

- ✓ Thermo Scientific™ ART™ Barrier Low-Retention Pipette Tips are specifically designed for eliminating sample retention within the pipette tip.
- ✓ Inner surface is extremely hydrophobic, increasing sample reproducibility and sample accuracy
- ✓ Crystal clear tips allow you to easily see transparent samples
- ✓ Tips seal when exposed to potential contaminants
- ✓ Sterile
- ✓ Certified RNase, DNase, ATP, Endotoxin, and Pyrogen Free
- ✓ Length : 7.9 cm
- ✓ Volume : 1000 µL

Sales unit :

800 Tips
