

**Product Name:** **Cultubes™**

**Catalogue No.:** **T406-3**

This document replaces any previous version

## 1. Product Description:

- **14 ml Culture Tube with Cap:** Disposable, sterile 17 x 95 mm culture tube with graduations and white marking area

## 2. Packaging:

- **Case:** 8 trays of 125 units / 1000 units per case

## 3. Product Specifications:

- **Material:**
  - Tube: Polystyrene
  - Cap: Polyethylene
- **Dual-position snap cap:**
  - Sealed position for anaerobic cultures
  - Loose position for aerobic work
- **Certified Pyrogen Free**
- **Gamma radiation sterilized at a SAL of  $10^{-3}$ ; specified dose between 6.5 kGy and 13.5 kGy**
- **Temperature range for tube only: -80°C to +80°C**
- **Not autoclavable**
- **Centrifuge resistant up to 1 400 g**

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resins conform to FDA 21 CFR 177.1640 & CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastics are in conformity to REACH standards
- **LATEX:** Materials are Latex Free

## 5. Quality Assurance:

- **Clear, no presence of contamination in plastic**
- **Visual attributes**
- **Dimension and volume measurements**
- **Closure verification**

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging which could compromise the product sterility

## 8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*