

Product Name: **Scintillation vial**

Catalogue No.: **S207**

This document replaces any previous version

1. Product Description:

- **20 ml Scintillation vial:** Non sterile, 6.5 ml vial

2. Packaging:

- **Case:** 1000 unit per case

3. Product Specifications:

- Material:
 - Tube: High Density Polyethylene
 - Cap: Polypropylene
- Temperature range: -80°C to +100°C
- Non-autoclavable

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions.

OWNERSHIP OF MATERIALS:

Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.