

Product Name: Dual Position Cap**Catalogue No.: T401-10DSPE**

This document replaces any previous version

1. Product Description:

- **Two-Position Snap Cap for 17 mm Polystyrene Culture Tube:** Non-sterile, disposable, and dual-position snap cap.
- **Available Color:**
 - Natural
- **Compatible with:**
 - T400-10

2. Packaging:

- **Case:** 1 bag of 1 000 units per case

3. Product Specifications:

- Material: Polyethylene
- Dual-position snap cap
 - Sealed position for anaerobic cultures
 - Loose position for aerobic work
- Non autoclavable

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG:** Plastics and colorants are in conformity with CONEG standards for heavy metals

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attribute
- Dimension measurements
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging.

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*