

**Product Name: DrainRack™ Jr.**

**Catalogue No.: M710-50 series**

This document replaces any previous version

## 1. Product Description:

- **Drain rack:** Slides are stored in an upright position for easy insertion and removal. Rack has 50 numbered slots that fits 50 to 100 slides.
- **Available Colors:**
  - Blue, Green, Pink White and Yellow
- **Compatible with:**
  - Fits standard 3 x 1 inch microscope slides

## 2. Packaging:

- **Case:** 10 individually wrapped units per case

## 3. Product Specifications:

- **Material:** High Impact Polystyrene
- **Temperature range:** -80°C to +80°C
- **Not Autoclavable**

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1640
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Closure verification

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

➤ **The lot number can be found in one or all of these locations:**

1. On exterior case label
2. On label inserted inside the master case
3. On the inner bag

**7. Storage Conditions:**

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging or product

**OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*