

Product Name: Histosette® II

Catalogue No.: M492-X Series

This document replaces any previous version

1. Product Description:

- **Tissue Cassette with Assembled Lid:** Non-Sterile, disposable cassette
- **Available Colors:**
 - White, Pink, Green, Yellow, Blue, Peach, Tan, Gray, Lilac, Orange and Aqua
 - Fluorescent Pink, Green and Yellow also available, minimum order quantity required

2. Packaging:

- **Case:** 3 inner packs of 500 units / 1 500 units per case

3. Product Specifications:

- **Material:**
 - Acetal
- **Writing surface angle:** 45°

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.2470
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free
- **CE:** Product is CE marked

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Lid closure verification

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper lid closure of the cassette before using
- Avoid picric acid, prolonged use in this acid will cause degradation
- Follow chemical resistance chart recommendations
- For use in automated equipment, follow the equipment manufacturer's instructions

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*