

**Product Name: OneHand™ Rack**

**Catalogue No.: T360-50**

This document replaces any previous version

## 1. Product Description:

- **Rack:** Non-sterile, 50 position, stackable rack with serrated holes for tubes with ribs and locking base for one-hand operations with skirted bottom Simport Microwtube® 0.5 and 2 ml
- **Available colors:**
  - Blue, Green, Lilac, Orange and Yellow
- **Compatible with:**
  - 0.5 to 2 ml Microwtubes®

## 2. Packaging:

- **Case:** 10 units per case

## 3. Product Specifications:

- **Material:**
  - Acetal
- Temperature range: -80°C to +30°C
- Not Autoclavable
- Total capacity: 50 Microtubes

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes

## 6. Traceability:

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;

## 8. Recommended Use:

Rack for holding tubes

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*