

Product Name: **Microcentrifuge Tube**

Catalogue No.: **T330-15 Series**

This document replaces any previous version

1. Product Description:

1.5ml Microcentrifuge Tube with Locking Cap: Non sterile, conical bottom, with attached Secure-Lock™ Pierceable cap

➤ **Available colors:**

- Natural, Blue, Green and Yellow

2. Packaging:

➤ **Case:** 10 packages of 500 units / 5 000 units per case
(packaged in tamperproof resealable safety-lock bags)

3. Product Specifications:

- **Material:**
 - Polypropylene
- Cap has a locking system to prevent opening during boiling applications
- Cap can be pierced easily with a syringe needle
- Certified RNase, DNase, Pyrogen and DNA Free
- Temperature range: -80°C to +121°C.
- Autoclavable at 121°C, for up to 30 minutes
- Molded graduation markings at 0.5ml, 1.0ml, and 1.5ml
- Etched surface on side of tube for sample identification
- Centrifuge resistant at up to 20 000g

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG / RoHS:** Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free
- **BSE / TSE:** Material is BSE / TSE Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Closure verification

6. Traceability:

- **Lot No. Composition:** 8 to 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- Verify proper cap closure when using biohazard material and / or chemical reagents;
- Follow chemical resistance chart recommendations;

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*