



# Simport

Since 1975

Scientific inc.

*A family owned Canadian company*

## SnapTwist™ Microwtube®

- **SAVE TIME** with a push-on twist-open microcentrifuge tube
- **NO MORE SCREW CAPS** to screw on and off
- **JUST PUSH** on the cap and voilà, tube is securely sealed
- **TWO TYPES OF CAPS** available



This revolutionary cap design is as tight and efficient as a screw cap with o-ring thanks to a deep internal lip. Readily usable for all routine lab applications including vortexing, centrifugation, freezing and gentle heating\*.

Push with your thumb to close, and turn just ¼ turn to open. No tightening required.

Caps are available separately in many colors as well as with an attachment loop.

**Don't take our word for it, try it for yourself.**  
Samples available on request.



**PUSH CAP DOWN  
FOR A SECURE SEAL**



**¼ TURN  
FAST REMOVAL**



**T342-XT**

Offered plain, skirted or conical bottom and non sterile.



**T343XLS  
Caps without loop**

**T343XLSL  
Caps with loop**

\* Not suitable for extended hard boiling.